

LIFE SCIENCE

# DP75

Digital Microscope Camera

Empowering Microscopic Discoveries



**EVIDENT**


Not for clinical diagnosis use.

# One Camera. Multiple Applications.


The DP75 digital microscope camera is a high-performance, multi-application imaging tool that makes it easy to capture high-resolution brightfield or fluorescence images using a single camera. It simplifies your microscopy imaging, so you can focus more on your work.




Colon. Stain: HE




Pig liver. Stain: Masson trichrome




Rat brain. Stain: DAPI Alexa Fluor 555 Cy5




Cytology. Stain: ThinPrep



Rat Pancreas. Stain: DAPI AF555 Cy5



Rat marrow. Stain: Giemsa



Mammary Gland. Stain: HER2

# See More without Switching Cameras

With the DP75 digital microscope camera, there is no need for switching between color and monochrome cameras when you want to capture brightfield and fluorescence images. Now, you can capture both in outstanding quality using a single powerful camera.

Designed to excel in a wide variety of applications, the DP75 camera offers a high-sensitivity cooled CMOS sensor, fluorescence images up to Cy 7.5 via a switchable infrared (IR) cut filter, and high-resolution imaging, making it a versatile tool that is up to the challenge of modern imaging tasks.


## Sharper Images, Clearer Insights

The DP75 camera makes it easier than ever to capture sharp, low-noise images. To push the image quality even further we integrated our real-time TruAI denoising algorithm into the camera. In addition, the DP75 maintains the high color fidelity our cameras are known for, delivering exceptional color reproduction and making your images as vivid as looking through the microscope oculars.

When imaging live specimens, a fast frame rate is important for efficiency and capturing the dynamics of your samples. With a fast frame rate of 22 frames per second (fps) at over 4K resolution and 60 fps at full HD resolution, the camera provides smooth, fast, live images for easy framing and comfortable live observation.




Mouse




BPAE cells

## High-Resolution, Wide Field of View Imaging

The camera's wide field of view capabilities enable you to find your target areas quickly, making your research more efficient. In addition, the DP75 camera enables you to capture high-resolution images even at low magnification with a maximum resolution of 8192 × 6000 through pixel shifting modes.




Without pixel shift



With pixel shift

# Quantitative Image Data with Wide-Wavelength Fluorescence


If you are using reagents with infrared emission for multicolor fluorescence imaging, the DP75 camera will transform the way you capture images.



Rat brain. Stain: DAPI Alexa Fluor 555 Cy5

## Infrared Fluorescence Imaging


The DP75 camera supports multiple staining combinations and wavelengths up to 1000 nm with a switchable IR cut filter. This setup enables you to—for example—check your sample conditions with your standard widefield microscope before spending time on a confocal microscope to finalize your imaging.



Ptk2. Stain: Cy7

## Qualitative Analysis Capabilities

The camera's linear mode enables intensity analysis without needing a dedicated monochrome camera. You also can access raw RGB pixel values for quantitative data regarding staining density or brightness. Moreover, the camera enables you to easily overlay fluorescence and brightfield images with pixel precision since you are using the same sensor for brightfield and fluorescence. This enables you to precisely identify the location of fluorescent expression, helping you focus on the relevant morphology and localization of your specimen.



Original: brightfield and fluorescence images of a zebrafish




Overlaid image

# Intelligent Features, Stunning Results

The DP75 camera can make your microscope observations more comfortable and efficient with smart features and AI support.


## Smart Observation Detection

The AI-based scene detection feature automatically recognizes five observation methods (brightfield, fluorescence, phase contrast, differential interference contrast, and polarization), enabling anyone to obtain high-quality images with minimal training.



AI-scene detection OFF

AI-scene detection ON




AI-scene detection OFF

AI-scene detection ON

## Multiple Image Alignment (MIA) Capabilities

The instant multiple image alignment (MIA) function simplifies the creation of wide-area images by moving the XY stage manually without any motorized setup, and the integrated position navigator helps ensure that you always know your position on the sample during brightfield and fluorescence imaging.



## Flexible Upgrades

Since the DP75 camera uses USB 3.1 Gen2, it is compatible with most PCs for a simple, effective upgrade to your current system.



\*Please see the system PC requirements in the specifications section.

## DP75 Specifications

| Item  |                        | Specifications   |
|---|------------------------|--|
| Camera type                                 |                        | Single-chip color camera (pixel shifting)<br>Cooling system: Peltier device (active cooling)   |
| Imaging sensor size                         |                        | 1.1-inch, 12.37-megapixel color CMOS image sensor, global shutter  |
| Camera mount                                |                        | C-mount  |
| Effective image resolution                  |                        | 8192 × 6000 (pixel shifting), 4096 × 3000 (3CMOS mode), 4096 × 3000 (1 × 1), 3840 × 2160 (1 × 1) (by cropped), 2048 × 1500 (2 × 2), 2048 × 1500 (1 × 1), 1920 × 1080 (1 × 1), ROI  |
| Sensitivity                                 |                        | 1x/2x/4x/8x/16x/32x<br>(ISO 100/200/400/800/1600/3200 equivalent)  |
| A/D   |                        | 12-bit   |
| Metering modes                              | Mode                   | Auto / SFL-Auto / Manual   |
|   | Adjustment             | ±2.0 EV step: 1/3 EV   |
|   | Time                   | 28 μs–120 s  |
| Binning                                     |                        | 2 × 2  |
| White balance                               |                        | Auto/One-touch/Manual/Area designation   |
| Black balance                               |                        | Auto/One-touch/Manual/Area designation   |
| Live frame rate *1                          |                        | 4096 × 3000 (1 × 1): 22 fps, 2048 × 1500 (2 × 2): 22 fps, 2048 × 1500 (1 × 1): 44 fps, 1920 × 1080 (1 × 1): 60 fps   |
| Still image transfer time                   |                        | 8192 × 6000 (pixel shifting): approx. 3 s, 4096 × 3000 (3CMOS Mode): approx. 2 s, 4096 × 3000 (1 × 1): approx. 1.2 s, 2048 × 1500 (2 × 2): approx. 1.0 s, 2048 × 1500 (1 × 1): approx. 0.4 s, 1920 × 1080 (1 × 1): approx. 0.4 s |
| Monochrome mode                             |                        | Available (Standard/Custom )   |
| Color space                                 |                        | sRGB, AdobeRGB*2   |
| Linear mode                                 |                        | Available  |
| IR cut filter                               |                        | Switchable: In: 400 nm up to 650 nm<br>Out: 400 nm up to 1000 nm   |
| Manual panoramic imaging (instant MIA) *3*4 |                        | Available (supports fluorescence as well as brightfield)   |
| Auto scene recognition mode*4               |                        | Available using an AI algorithm (supports: brightfield, fluorescence, phase contrast, differential interference contrast, and polarization)  |
| Position navigator*4                        |                        | Available  |
| Control software                            |                        | cellSens Entry / Standard / Dimension v. 4.2.1 or later<br>DP2-TWAIN v. 10.5. or later   |
| External trigger                            |                        | Available (input/output)   |
| Dimensions (W × D × H)                      | Camera interface cable | Approx. 2.7 m (8.9 ft)   |
|   | AC adapter             | 107 × 47 × 30 mm (4.2 × 1.9 × 1.2 in.)/Approx 0.3 kg (0.7 lb)  |

## DP75 System Requirements

| Item   | Specifications   |
|--------|--|
| CPU    | Intel Core i5, Intel Core i7, Intel Xeon, or equivalent of Intel CPU |
| RAM    | 8 GB or more (recommended 16 GB or more)                             |
| PC I/F | USB 3.1 Gen2 (TypeA) (a dedicated board is unnecessary)*5            |
| OS     | Windows 10 Pro (64-bit)  |
|        | Windows 11 Pro (64-bit)  |

\*1 Frame rate may decrease depending on the condition of your PC, monitor resolution, and/or software.


\*2 Monitor designed to meet Adobe RGB is required.

\*3 Manual Process option license is required for cellSens Standard.


\*4 Not available in the combination of cellSens Entry or DP2-TWAIN.

\*5 Operatable with USB3.1 Gen1 (5 Gbps), but framerate is decreased.

## System Diagram



## DP75 Dimensions



Cover image: Rat brain. DAPI, Alexa Fluor 488, Alexa Fluor 555, Alexa Fluor 750

- **EVIDENT CORPORATION is ISO14001 certified.**  
For details on certification registration, visit <https://www.olympus-lifescience.com/en/support/iso/>
- **EVIDENT CORPORATION is ISO9001 certified.**
- All company and product names are registered trademarks and/or trademarks of their respective owners.
- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.
- Images on the PC monitors are simulated.

[EvidentScientific.com](https://www.evidentscientific.com)

**EVIDENT**

EVIDENT CORPORATION  
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0910, Japan

**OLYMPUS**