


Virtual Slide Management




Net Image Server^{SQL}*

The Net Image Server^{SQL} for the Olympus VS120®* system is a versatile image access and storage solution for digital virtual microscopy. This client-server-based data management system ensures that even with large amounts of data and multiple users, workflows run efficiently with outstanding data security. Since storage and traceability of images, reports and analytical results is vital, the platform incorporates a clear, easy-to-manage structure. The image data remains on the system, while only the required image is downloaded to the client computer; ensuring smooth data access is maintained, even when vast numbers of users are online. Images can easily be accessed, viewed and annotated with an internet browser or using OlyVIA®, free-of-charge image viewer. Each user can review the sample slide and control the position and magnifications, as if were a real slide. The Net Image Server^{SQL} comes with integrated security settings, giving tutors controlled access to image collections and related data.

Net Image Server^{SQL}

Features and Benefits

- Stable, secure and dependable
- Outstanding speed of virtual slide image access
- Supports Microsoft® SQL Server 2003 (R2 32Bit), MS SQL 2005 Express/Standard and MS SQL 2008 Express/Standard
- Password protection for access and management of digital
- Standard authentication via SQL or Windows® authentication login
- Supports customized database fields and user defined database structure
- Functionality beyond virtual slides; Microsoft Word, Excel, Powerpoint documents as well as macro images captured by other devices
- Provides web access for images and conferencing capabilities
- Supports a wide range of server hardware configurations



Slide Into Digital™ - Powered by Olympus Net Image Server^{SQL}

Olympus has always provided products of exceptional quality and outstanding performance in the medical and veterinary educational fields. From our educational microscopes, to the VS120 Virtual Slide System, Olympus has been at the forefront of providing innovative tools for today's top educators and their students. With Slide into Digital, our new educational platform powered by Olympus Net Image Server^{SQL}, we continue this legacy by providing data management and image sharing solutions for today's digital classroom.

With Slide Into Digital

- Create a standardized, uniform collection of virtual slides that will never deteriorate or break
- Improve your curriculum and excite students with innovative, popular technology
- Free OlyVIA viewing software for all your students
- Students can access slides 24 hours a day, 7 days a week on or off campus
- No custom design work needed
- Reliable and durable server software insuring up time, even during peak activity periods
- Annual or One-Time licenses available
Remote Installation and training included
- Convert your old Bacus webslides to Olympus .vsi with a 30-day license to VS-Desktop®
- Over 3000 virtual slides available for your educational program in the Olympus Virtual Slide Library

VS Desktop, Slide Into Digital, VS and OlyVIA is a trademark of Olympus America Inc.
All other third party marks are property of their respective owners.

VS120, OlyVIA and Net Image Server are not for clinical diagnostic use in the United States.

OLYMPUS®

OLYMPUS AMERICA INC.
3500 Corporate Parkway, Center Valley, PA 18034-0610, U.S.A.

OLYMPUS CANADA INC.
25 Leek Crescent, Richmond Hill, Ontario L4B 4B3

OLYMPUS LATIN AMERICA INC.
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126, U.S.A.